

Programming data for Spectroquant® test kits for photometer SQ 118

The programming data until program place 174 correspond with version 1.80 of the photometer SQ 118. The other methods are not stored in the memory of the photometer SQ 118. These methods can be stored manually in the memory of the photometer SQ 118.

Curve = SQ 118 contains the correct pre-programmed calibration curve as a non-linear function

CT = Cell Test

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor	
000	Absorption		0.000	-	3.000		A					
001	Aluminium	Al	14825	0.05	-	0.35	mg/l	Al	2+0	50	550	Curve
002	Aluminium	Al	14825	0.20	-	1.50	mg/l	Al	2+0	10	550	0.550
003	Aluminium	Al	14825	2.0	-	13.0	mmol/m ³	Al	2+0	50	550	Curve
004	Aluminium	Al	14825	7.4	-	55.6	mmol/m ³	Al	2+0	10	550	20.4
005	Ammonium	NH ₄	14752	0.010	-	0.800	mg/l	NH ₄	5+5	50	690	0.317
006	Ammonium	NH ₄	14752	0.10	-	3.50	mg/l	NH ₄	5+5	10	690	1.58
007	Ammonium	NH ₄ -N	14752	0.010	-	0.620	mg/l	NH ₄ -N	5+5	50	690	0.246
008	Ammonium	NH ₄ -N	14752	0.10	-	2.70	mg/l	NH ₄ -N	5+5	10	690	1.23
009	Ammonium	NH ₄	14752	0.6	-	44.4	mmol/m ³	NH ₄	5+5	50	690	17.6
010	Ammonium	NH ₄	14752	6	-	195	mmol/m ³	NH ₄	5+5	10	690	87.6
011	CT Ammonium	NH ₄	14544	Please use new programming data from page 7								
012	CT Ammonium	NH ₄ -N	14544									
013	CT Ammonium	NH ₄	14544									
014	Boron	B	14839	0.050	-	0.800	mg/l	B	12+2	10	550	0.323
015	Calcium	Ca	14815	5	-	160	mg/l	Ca	3+0	20	550	68.5
016	Calcium	CaO	14815	7	-	224	mg/l	CaO	3+0	20	550	95.8
017	Calcium	CaCO ₃	14815	0.13	-	4.00	mol/m ³	CaCO ₃	3+0	20	550	1.71
018	Chlorine	Cl ₂	14828	Item deleted, replacements see chlorine tests at page 9								
019	Chlorine	Cl ₂	14828									
020	Chloride	Cl	14755	1.0	-	20.0	mg/l	Cl	0+0	20	446	9.23
021	Chloride	Cl	14755	10	-	150	mg/l	Cl	0+0	10	525	Curve
022	Chloride	Cl	14755	28	-	564	mmol/m ³	Cl	0+0	20	446	260
023	Chloride	Cl	14755	0.28	-	4.23	mol/m ³	Cl	0+0	10	525	Curve

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor	
024	Chromium	Cr	14758	0.010	-	0.600	mg/l	Cr	0+0	50	550	0.254
025	Chromium	Cr	14758	0.10	-	3.00	mg/l	Cr	0+0	10	550	1.27
026	Chromium	Cr	14758	0.2	-	11.5	mmol/m ³	Cr	0+0	50	550	4.88
027	Chromium	Cr	14758	2.0	-	58.0	mmol/m ³	Cr	0+0	10	550	24.4
028	CT COD	COD	14540	Please use new programming data from page 10								
029	CT COD	COD	14541	Please use new programming data from page 10								
030	Cyanide	CN	14800	0.002	-	0.100	mg/l	CN	5+0	50	585	0.0416
031	Cyanide	CN	14800	0.025	-	0.500	mg/l	CN	5+0	10	585	0.208
032	Iron	Fe	14761	0.01	-	1.00	mg/l	Fe	3+0	50	565	0.420
033	Iron	Fe	14761	0.20	-	2.50	mg/l	Fe	3+0	20	565	1.05
034	Iron	Fe	14761	0.2	-	18.0	mmol/m ³	Fe	3+0	50	565	7.52
035	Iron	Fe	14761	3.6	-	45.0	mmol/m ³	Fe	3+0	20	565	18.8
036	Gold	Au	14821	0.5	-	12.0	mg/l	Au	0+0	10	550	6.06
037	Hydrazine	N ₂ H ₄	14797	Item deleted, replacement is item 09711, see page 10								
038	Hydrazine	N ₂ H ₄	14797									
039	Hydrazine	N ₂ H ₄	14797									
040	Hydrazine	N ₂ H ₄	14797									
041	Copper	Cu	14767	0.05	-	2.00	mg/l	Cu	5+0	50	585	0.910
042	Copper	Cu	14767	0.20	-	6.00	mg/l	Cu	5+0	10	585	4.55
043	Copper	Cu	14767	0.8	-	31.5	mmol/m ³	Cu	5+0	50	585	14.3
044	Copper	Cu	14767	3	-	126	mmol/m ³	Cu	5+0	10	585	71.6
045	Manganese	Mn	14770	0.01	-	2.50	mg/l	Mn	2+2	50	446	1.12
046	Manganese	Mn	14770	0.5	-	10.0	mg/l	Mn	2+2	10	446	5.62
047	Manganese	Mn	14770	0.2	-	46.0	mmol/m ³	Mn	2+2	50	446	20.5
048	Manganese	Mn	14770	9	-	182	mmol/m ³	Mn	2+2	10	446	102
049	Nickel	Ni	14785	0.02	-	2.00	mg/l	Ni	1+2	50	446	0.952
050	Nickel	Ni	14785	0.20	-	5.00	mg/l	Ni	1+2	20	446	2.38
051	Nickel	Ni	14785	0.4	-	34.0	mmol/m ³	Ni	1+2	50	446	16.2
052	Nickel	Ni	14785	3.0	-	85.0	mmol/m ³	Ni	1+2	20	446	40.6



Programming data for Spectroquant® test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor	
053	Nitrate	NO ₃	14773	1.0	-	25.0	mg/l	NO ₃	10+0	50	525	8.35
054	Nitrate	NO ₃	14773	5.0	-	90.0	mg/l	NO ₃	10+0	10	525	41.7
055	Nitrate	NO ₃ -N	14773	0.20	-	5.60	mg/l	NO ₃ -N	10+0	50	525	1.89
056	Nitrate	NO ₃ -N	14773	1.1	-	20.3	mg/l	NO ₃ -N	10+0	10	525	9.43
057	Nitrate	NO ₃	14773	16	-	403	mmol/m ³	NO ₃	10+0	50	525	135
058	Nitrate	NO ₃	14773	0.08	-	1.45	mol/m ³	NO ₃	10+0	10	525	0.672
059	CT Nitrate	NO ₃	14542	2.0	-	80.0	mg/l	NO ₃	10+0	16	525	32.8
060	CT Nitrate	NO ₃ -N	14542	0.5	-	18.0	mg/l	NO ₃ -N	10+0	16	525	7.41
061	CT Nitrate	NO ₃	14542	32	-	1290	mmol/m ³	NO ₃	10+0	16	525	529
062	Nitrite	NO ₂	14776	Please use new programming data from page 13								
063	Nitrite	NO ₂	14776									
064	Nitrite	NO ₂ -N	14776									
065	Nitrite	NO ₂ -N	14776									
066	Nitrite	NO ₂	14776									
067	Nitrite	NO ₂	14776									
068	CT Nitrite	NO ₂	14547	0.05	-	2.00	mg/l	NO ₂	10+0	16	525	0.920
069	CT Nitrite	NO ₂ -N	14547	0.020	-	0.610	mg/l	NO ₂ -N	10+0	16	525	0.280
070	CT Nitrite	NO ₂	14547	1.1	-	43.5	mmol/m ³	NO ₂	10+0	16	525	20.0
071	P(PMB)	P	14848	Please use new programming data from page 14								
072	P(PMB)	P	14848									
073	P(PMB)	PO ₄	14848									
074	P(PMB)	PO ₄	14848									
075	P(PMB)	P ₂ O ₅	14848									
076	P(PMB)	P ₂ O ₅	14848									
077	P(PMB)	PO ₄	14848									
078	P(PMB)	PO ₄	14848									
079	CT P(PMB)	P	14543	0.05	-	5.00	mg/l	P	5+0	16	690	1.64
080	CT P(PMB)	PO ₄	14543	0.2	-	15.3	mg/l	PO ₄	5+0	16	690	5.03
081	CT P(PMB)	P ₂ O ₅	14543	0.1	-	11.5	mg/l	P ₂ O ₅	5+0	16	690	3.76

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor	
082	CT P(PMB)	PO ₄	14543	2	-	160	mmol/m ³	PO ₄	5+0	16	690	53.0
083	P (VM)	P	14842	1.0	-	30.0	mg/l	P	0+0	10	405	16.1
084	P (VM)	PO ₄	14842	3.0	-	92.0	mg/l	PO ₄	0+0	10	405	49.4
085	P (VM)	P ₂ O ₅	14842	2.0	-	69.0	mg/l	P ₂ O ₅	0+0	10	405	36.9
086	P (VM)	PO ₄	14842	32	-	970	mmol/m ³	PO ₄	0+0	10	405	520
087	CT P (VM)	P	14546	0.5	-	25.0	mg/l	P	0+0	16	405	12.2
088	CT P (VM)	PO ₄	14546	1.5	-	77.0	mg/l	PO ₄	0+0	16	405	37.4
089	CT P (VM)	P ₂ O ₅	14546	1.1	-	57.0	mg/l	P ₂ O ₅	0+0	16	405	28.0
090	CT P (VM)	PO ₄	14546	16	-	807	mmol/m ³	PO ₄	0+0	16	405	394
091	Silicate	Si	14794	0.010	-	0.800	mg/l	Si	3+5	50	820	0.255
092	Silicate	Si	14794	0.10	-	5.00	mg/l	Si	3+5	20	660	1.79
093	Silicate	SiO ₂	14794	0.02	-	1.70	mg/l	SiO ₂	3+5	50	820	0.545
094	Silicate	SiO ₂	14794	0.2	-	10.7	mg/l	SiO ₂	3+5	20	660	3.83
095	Silicate	SiO ₂	14794	0.3	-	28.5	mmol/m ³	SiO ₂	3+5	50	820	9.08
096	Silicate	SiO ₂	14794	3	-	178	mmol/m ³	SiO ₂	3+5	20	660	63.7
097	CT Sulfate	SO ₄	14548	Please use new programming data from page 16								
098	CT Sulfate	SO ₄	14548									
099	Sulfide	S	14779	0.020	-	0.500	mg/l	S	1+0	50	660	0.260
100	Sulfide	S	14779	0.25	-	4.00	mg/l	S	1+0	10	660	Curve
101	Sulfide	S	14779	0.6	-	15.1	mmol/m ³	S	1+0	50	660	7.86
102	Sulfide	S	14779	7	-	121	mmol/m ³	S	1+0	10	660	Curve
103	Zinc	Zn	14832	Please use new programming data from page 18								
104	Zinc	Zn	14832									
105	CT Zinc	Zn	14566	0.20	-	5.00	mg/l	Zn	15+0	16	495	4.88
106	CT Chromium	Cr	14552	0.05	-	2.00	mg/l	Cr	0+0	16	550	0.971
107	CT Iron	Fe	14549	0.10	-	4.00	mg/l	Fe	3+0	16	565	1.63
108	CT Copper	Cu	14553	Please use new programming data from page 10								
109	CT Nickel	Ni	14554	0.10	-	6.00	mg/l	Ni	1+2	16	446	3.85
110	CT Cyanide	CN	14561	0.025	-	0.500	mg/l	CN	10+0	16	585	0.313

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor
111	CT Phenol	PHE 14551	0.10	-	2.50	mg/l	PHE	1+0	16	495	3.39
112	CT COD	COD 14555	500	-	10000	mg/l	COD	0+0	16	585	4630
113	Turbidity		10	-	400	TE(F)		0+0	50	525	145
114	CT Lead	Pb 14833	Please use new programming data from page 11								
115	CT Cadmium	Cd 14834	Please use new programming data from page 8								
116	CT Fluoride	F 14557	Please use new programming data from page 10								
117	CT Nitrate	NO ₃ 14556	0.5	-	15.0	mg/l	NO ₃	15+0	16	495	6.50
118	CT Nitrate	NO ₃ -N 14556	0.11	-	3.40	mg/l	NO ₃ -N	15+0	16	495	1.47
119	CT Nitrate	NO ₃ 14556	8	-	242	mmol/m ³	NO ₃	15+0	16	495	105
120	CT HCHO	HCHO 14500	0.1	-	10.0	mg/l	HCHO	5+0	16	585	3.36
121	Sulfate	SO ₄ 14791	25	-	300	mg/l	SO ₄	5+7	20	525	122
122	Iron	Fe 8023	0.2	-	20.0	mg/l	Fe	10+0	10	525	7.00
123	Silver	Ag 14831	Please use new programming data from page 15								
124	Silver	Ag 14831									
125	Urea	Urea 14544	0.5	-	25.0	mg/l		5+5	16	690	18.0
126	CT Ammonium	NH ₄ 14558	0.2	-	10.0	mg/l	NH ₄	15+0	16	690	5.71
127	CT Ammonium	NH ₄ -N 14558	0.20	-	7.80	mg/l	NH ₄ -N	15+0	16	690	4.43
128	CT Ammonium	NH ₄ 14559	5	-	100	mg/l	NH ₄	15+0	16	690	51.1
129	CT Ammonium	NH ₄ -N 14559	4.0	-	78.0	mg/l	NH ₄ -N	15+0	16	690	39.7
130	CT Nitrate	NO ₃ 14563	2	-	100	mg/l	NO ₃	10+0	16	340	65.1
131	CT Nitrate	NO ₃ -N 14563	0.5	-	22.6	mg/l	NO ₃ -N	10+0	16	340	14.7
132	CT COD	COD 14560	4.0	-	40.0	mg/l	COD	0+0	16	340	-41.0
133	CT Sulfate	SO ₄ 14564	Please use new programming data from page 16								
134	CT Sulfite	SO ₃ 14394	1.0	-	25.0	mg/l	SO ₃	2+0	16	405	9.09
135	Iodine colour no	IFZ	0.05	-	3.00	IFZ		0+0	10	340	1.27
136	Iodine colour no	IFZ	1.0	-	50.0	IFZ		0+0	10	446	17.0
137	Hazen colour no	Hazen	2	-	200	Hazen		0+0	50	340	78.0
138	Hazen colour no	Hazen	10	-	1000	Hazen		0+0	10	340	377
139	Colour	1/m	0.5	-	50.0	1/m		0+0	50	446	20.0

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor	
140	CT Total	N	14537	0.5	-	18.0	mg/l	N	10+0	16	525	8.00
141	CT Potassium	K	14562	Please use new programming data from page 15								
142	CT Phosphorus	P	14729	1.0	-	25.0	mg/l	P	5+0	16	690	8.33
143	CT Phosphorus	PO ₄	14729	3.1	-	76.8	mg/l	PO ₄	5+0	16	690	25.5
144	CT Surfactants anionic	Surf. anion.	14697	0.05	-	2.00	mg/l	MBAS	5+0	16	660	2.27
145	CT Hydrogen peroxide	H ₂ O ₂	14731	2.0	-	20.0	mg/l	H ₂ O ₂	2+0	16	405	40.0
146	Chlorine dioxide	ClO ₂	14828	Item deleted, replacement is item 00608, see page 9								
147	Ozone	O ₃	14828	Item deleted, replacement is item 00607, see page 13								
148	CT C OD	COD	14690	50	-	500	mg/l	COD	0+0	16	446	-400
149	CT COD	COD	14691	300	-	3500	mg/l	COD	0+0	16	585	3280
150	CT Chloride	Cl	14730	Please use new programming data from page 9								
151	CT Magnesium	Mg	14684	Item deleted, replacement is item 00815, see page 11								
152	CT Tin	Sn	14622	0.10	-	2.50	mg/l	Sn	15+0	16	660	1.22
153	CT Total hardness	TH	14565	Item deleted, replacement is item 00961, see page 17+18								
154	CT Hardness Ca	Ca	14565									
156	CT Residual hardness	RH	14683	0.075	-	0.750	°d		10+0	16	565	0.553
157	CT Residual hardness	Ca	14683	0.50	-	5.00	mg/l	Ca	10+0	16	565	3.95
158	HCHO	HCHO	14678	0.02	-	1.80	mg/l	HCHO	10+0	50	585	0.82
159	HCHO	HCHO	14678	1.00	-	9.00	mg/l	HCHO	10+0	10	585	4.10
160	Palladium	Pd		0.10	-	1.50	mg/l	Pd	1+0	10	525	0.57
161	Platinum	Pt		0.10	-	1.50	mg/l	Pt	5+1	10	690	2.36
162	Mercury	Hg		0.025	-	1.000	mg/l	Hg	5+0	50	565	0.563
163	Antimony	Sb		0.25	-	8.00	mg/l	Sb	3+2	10	635	4.44
164	Arsenic	As		Application deleted, replacement is item 01747, see page 8								

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor	
165	Chromium Bath	CrO ₃	20	-	400	g/l	CrO ₃	0+0	10	446	540	
166	Copper Bath	Cu	10.0	-	80.0	g/l	Cu	0+0	10	820	42.0	
167	Nickel Bath	Ni	10	-	120	g/l	Ni	0+0	10	690	60.0	
168	CT Oxygen	O ₂	14694	0.5	-	12.0	mg/l	O ₂	0+0	16	495	12.8
169	Chlorine dioxide	ClO ₂	14732	Item deleted, replacement is item 00608, see page 9								
170	Chlorine dioxide	ClO ₂	14732									
171	Chlorine	Cl ₂	14732									
172	Chlorine	Cl ₂	14732	Item deleted, replacements see chlorine tests at page 9								
173	Ozone	O ₃	14732									
174	Ozone	O ₃	14732	Item deleted, replacement is item 00607, see page 13								
CT Acid capacity to pH 4.3	OH	01762	Item deleted									
CT Acid capacity to pH 4.3		01758	can not be measured in the SQ 118 because of not available wavelength									
CT Alcohol	C ₂ H ₅ OH	14965	Test deleted. Reflectoquant Test is replacement. Additionally RQflex required.									
Alkalinity	see Acid capacity to pH 4.3											
CT Aluminium	Al	00594	0.02	-	0.50	mg/l	Al	5+0	16	550	non linear	
										factor a	0	
										factor b	0.8021	
										factor c	-0.6126	
										factor d	0.3009	
										factor e	0	
CT Ammonium	NH ₄	14739	0.01	-	2.00	mg/l	NH ₄	15+0	16	690	1.13	
CT Ammonium	NH ₄ -N	14739	0.01	-	2.60	mg/l	NH ₄ -N	15+0	16	690	0.877	
CT Ammonium	NH ₄	14544	0.6	-	20.6	mg/l	NH ₄	15+0	16	690	10.1	
CT Ammonium	NH ₄ -N	14544	0.5	-	16.0	mg/l	NH ₄ -N	15+0	16	690	7.84	

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor
Ammonium	NH ₄ -N	00683	2.0	-	75.0	mg/l	NH ₄ -N	15+0	10	690	27.8
Ammonium	NH ₄ -N	00683	5	-	150	mg/l	NH ₄ -N	15+0	10	690	55.6
Ammonium	NH ₄	00683	2.6	-	96.4	mg/l	NH ₄	15+0	10	690	35.9
Ammonium	NH ₄	00683	6	-	193	mg/l	NH ₄	15+0	10	690	71.8
CT AOX	AOX	00675	0.05	-	2.50	mg/l	AOX	15+0	16	446	4.0
Arsenic	As	01747	0.001	-	0.020	mg/l	As	120+0	20	525	0.057
Arsenic	As	01747	0.005	-	0.100	mg/l	As	120+0	10	525	0.115
CT BOD	BOD	00687	0.5	-	12.0	mg/l	BOD	0+0	16	495	12.8
CT Boron	B	00826	0.05	-	2.00	mg/l	B	60+0	16	405	1.41
Bromine	Br ₂	00605	0.020	-	2.000	mg/l	Br ₂	1+0	50	550	1.50
Bromine	Br ₂	00605	0.10	-	10.00	mg/l	Br ₂	1+0	10	550	7.52
CT Cadmium	Cd	14834	0.025	-	1.000	mg/l	Cd	2+0	16	525	0.630
Cadmium	Cd	01745	0.002	-	0.100	mg/l	Cd	2+0	50	525	0.141
Cadmium	Cd	01745	0.010	-	0.500	mg/l	Cd	2+0	10	525	0.704
Calcium ¹⁾	Ca	00049	0.20	-	4.00	mg/l	Ca	5+0	10	565	
CT Calcium	Ca	00858	10	-	250	mg/l	Ca	3+0	16	565	non linear
										factor a	0
										factor b	409.29
										factor c	-295.91
										factor d	138.96
										factor e	0

¹⁾ Remark to Calcium Test:

For this analysis an own calibration is required. The own calibration should be regularly check with a standard but at least every 3 month re-checked if everything is working fine. You will find all detailed information in our application list on the Internet homepage <http://photometry.merck.de>.



Programming data for Spectroquant® test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor
CT Calcium	CaCO3	00858	25	-	625	mg/l	CaCO3	3+0	16	565	non linear
										565	factor a
										565	factor b
										565	factor c
										565	factor d
										565	factor e
CT Calcium	CaO	00858	14	-	350	mg/l	CaO	3+0	16	565	non linear
										565	factor a
										565	factor b
										565	factor c
										565	factor d
										565	factor e
Carbohydrazide	see Oxygen Scavengers										
CT Chloride	Cl	14730	5	-	125	mg/l	Cl	0+0	16	525	113
Chloride	Cl	14897	2.5	-	25.0	mg/l	Cl	3+0	10	446	28.2
Chloride	Cl	14897	10	-	250	mg/l	Cl	3+0	10	495	107
CT Chlorine	Cl ₂	00595	0.03	-	6.00	mg/l	Cl ₂	1+0	16	550	2.47
CT Chlorine	Cl ₂	00597	0.03	-	6.00	mg/l	Cl ₂	1+0	16	550	2.47
Chlorine	Cl ₂	00598	0.010	-	1.000	mg/l	Cl ₂	1+0	50	550	0.666
Chlorine	Cl ₂	00598	0.05	-	6.00	mg/l	Cl ₂	1+0	10	550	3.33
Chlorine	Cl ₂	00599	0.010	-	1.000	mg/l	Cl ₂	1+0	50	550	0.666
Chlorine	Cl ₂	00599	0.05	-	6.00	mg/l	Cl ₂	1+0	10	550	3.33
Chlorine	Cl ₂	00602	0.010	-	1.000	mg/l	Cl ₂	1+0	50	550	0.666
Chlorine	Cl ₂	00602	0.05	-	6.00	mg/l	Cl ₂	1+0	10	550	3.66
Chlorine dioxide	ClO ₂	00608	0.020	-	2.000	mg/l	ClO ₂	2+1	50	550	1.27
Chlorine dioxide	ClO ₂	00608	0.10	-	10.00	mg/l	ClO ₂	2+1	10	550	6.33

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text	Item	Measuring range					Citation form	Reaction time	cells. mm	filter nm	factor
CT COD	COD	14540	10	-	150	mg/l	COD	0+0	16	446	-201
CT COD	COD	14541	25	-	1500	mg/l	COD	0+0	16	585	1670
CT COD	COD	14895	15	-	300	mg/l	COD	0+0	16	446	-222
CT COD (Hg free)	COD	09772	10	-	150	mg/l	COD	0+0	16	446	-238
CT COD (Hg free)	COD	09773	100	-	1500	mg/l	COD	0+0	16	585	1925
CT Copper	Cu	14553	0.05	-	8.00	mg/l	Cu	5+0	16	585	3.60
Cyanuric Acid ¹⁾	CYS	19250	2	-	160	mg/l	CYS	0+0	20	525	non linear
										factor a	-0,0757
										factor b	80,335
										factor c	3,7186
										factor d	0
										factor e	0
DEHA Diethylhydroxylamine	see Oxygen Scavengers										
Detergents	see Surfactants										
Dissolved Oxygen	see Oxygen										
KT Fluoride	F	14557	0.10	-	1.50	mg/l	F	5+0	16	635	1.73
Fluoride	F	14598	0.10	-	2.00	mg/l	F	5+0	10	635	2.40
Fluoride	F	14598	1.0	-	20.0	mg/l	F	5+0	10	635	25.7
Hardness	see Residual hardness or Total hardness										
Hydrazine	N ₂ H ₄	09711	0.005	-	0.400	mg/l	N ₂ H ₄	5+0	50	446	0.174
Hydrazine	N ₂ H ₄	09711	0.02	-	2.00	mg/l	N ₂ H ₄	5+0	10	446	0.870
Hydroquinone	see Oxygen Scavengers										
Hydrogen peroxide	H ₂ O ₂	18789	0.015	-	3.000	mg/l	H ₂ O ₂	10+0	20	446	1.33
Hydrogen peroxide	H ₂ O ₂	18789	0.03	-	6.00	mg/l	H ₂ O ₂	10+0	10	446	2.65

¹⁾ Remark to Cyanuric Acid Test:

As blank dist. water without adding reagents must be used.

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text	Item	Measuring range	Citation form	Reaction time	cells. mm	filter nm	factor				
Hydrogen sulfide	see Sulfide										
Iodine	I ₂	00606	0.050	-	2.000	mg/l	I ₂	1+0	50	550	2.38
Iodine	I ₂	00606	0.20	-	10.00	mg/l	I ₂	1+0	10	550	11.91
CT Iron	Fe	14896	1.0	-	50.0	mg/l	Fe	5+0	16	525	28.3
Iron	Fe	00796	0.010	-	1.000	mg/l	Fe	5+10	50	495	1.12
Iron	Fe	00796	0.10	-	5.00	mg/l	Fe	5+10	10	495	5.62
Isoascorbic acid (erythorbic acid)	see Oxygen Scavengers										
CT Lead	Pb	14833	0.10	-	5.00	mg/l	Pb	0+0	16	525	4.39
Lead	Pb	09717	0.010	-	1.000	mg/l	Pb	0+0	50	525	1.25
Lead	Pb	09717	0.10	-	5.00	mg/l	Pb	0+0	10	525	6.25
CT Magnesium	Mg	00815	5.0	-	75.0	mg/l	Mg	3+0	16	565	non linear
										factor a	0
										factor b	268.51
										factor c	-419.66
										factor d	475.72
										factor e	0
CT Manganese	Mn	00816	0.10	-	5.00	mg/l	Mn	2+5	16	446	4.13
Manganese	Mn	01739	0.005	-	0.400	mg/l	Mn	10+0	50	565	0.344
Manganese	Mn	01739	0.05	-	2.00	mg/l	Mn	10+0	10	565	1.72
Methylethylketoxime (2-butanoneoxime)	see Oxygen Scavengers										
CT Molybdenum	Mo	00860	can not be measured in the SQ 118 because of not available wavelength								
Molybdenum	Mo	19252	can not be measured in the SQ 118 because of not available wavelength								



Programming data for Spectroquant[®] test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor
Monochloramine	Cl ₂	01632	0.050	-	2.000	mg/l	Cl ₂	5+10	50	690	1.005
Monochloramine	Cl ₂	01632	0.25	-	10.00	mg/l	Cl ₂	5+10	10	690	5.025
Monochloramine	NH ₂ Cl	01632	0.036	-	1.450	mg/l	NH ₂ Cl	5+10	50	690	0.73
Monochloramine	NH ₂ Cl	01632	0.18	-	7.25	mg/l	NH ₂ Cl	5+10	10	690	3.65
Monochloramine	NH ₂ Cl-N	01632	0.010	-	0.392	mg/l	NH ₂ Cl-N	5+10	50	690	0.1985
Monochloramine	NH ₂ Cl-N	01632	0.05	-	1.96	mg/l	NH ₂ Cl-N	5+10	10	690	0.993
CT Nitrate	NO ₃ -N	14764	1.0	-	50.0	mg/l	NO ₃ -N	10+0	16	340	27.9
CT Nitrate	NO ₃	14764	4	-	220	mg/l	NO ₃	10+0	16	340	124
CT Nitrate	NO ₃ -N	00614	23	-	225	mg/l	NO ₃ -N	10+0	16	340	132.5
CT Nitrate	NO ₃	00614	100	-	1000	mg/l	NO ₃	10+0	16	340	586.7
Nitrate	NO ₃ -N	09713	0.10	-	5.00	mg/l	NO ₃ -N	10+0	50	340	3.92
Nitrate	NO ₃ -N	09713	1.0	-	25.0	mg/l	NO ₃ -N	10+0	10	340	19.6
Nitrate	NO ₃	09713	0.5	-	22.2	mg/l	NO ₃	10+0	50	340	17.4
Nitrate	NO ₃	09713	4.4	-	111.0	mg/l	NO ₃	10+0	10	340	86.8
Nitrate in seawater	NO ₃ -N	14942	0.2	-	17.0	mg/l	NO ₃ -N	15+60	10	495	non linear
										factor a	0
										factor b	9.73
										factor c	-1.83
										factor d	0.35
										factor e	0
Nitrate in seawater	NO ₃	14942	1.0	-	75.0	mg/l	NO ₃	15+60	10	495	non linear
										factor a	0
										factor b	43.38
										factor c	-8.53
										factor d	1.63
										factor e	0

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text	Item	Measuring range	Citation form	Reaction time	cells. mm	filter nm	factor				
CT Nitrit	NO ₂ -N	00609	1.0 - 90.0 mg/l	NO ₂ -N	20+0	16	585	78.7			
CT Nitrit	NO ₂	00609	3.3 - 295.2 mg/l	NO ₂	20+0	16	585	258.38			
Nitrite	NO ₂	14776	0.007 - 0.650 mg/l	NO ₂	10+0	50	525	0.248			
Nitrite	NO ₂	14776	0.07 - 3.25 mg/l	NO ₂	10+0	10	525	1.24			
Nitrite	NO ₂ -N	14776	0.002 - 0.200 mg/l	NO ₂ -N	10+0	50	525	0.0751			
Nitrite	NO ₂ -N	14776	0.02 - 1.00 mg/l	NO ₂ -N	10+0	10	525	0.376			
Nitrogen, Ammonium	see Ammonium										
Nitrogen, Nitrate	see Nitrate										
Nitrogen, Nitrite	see Nitrite										
CT Nitrogen, total	N	00613	0.5 - 15.0 mg/l	N	10+0	16	340	15.3			
CT Nitrogen, total	N	14763	10 - 150 mg/l	N	10+0	16	340	154			
Organic Acids, Volatile	see Volatile organic acids										
Organic Carbon,	see TOC										
Oxygen Scavengers	DEHA	19251	0.020 - 0.500 mg/l	DEHA	10+0	20	565	0.408			
Oxygen Scavengers	Carbohy	19251	0.027 - 0.667 mg/l	Carbohy	10+0	20	565	0.544			
Oxygen Scavengers	Hydro	19251	0.053 - 1.315 mg/l	Hydro	2+0	20	565	1.073			
Oxygen Scavengers	ISA	19251	0.078 - 1.950 mg/l	ISA	10+0	20	565	1.592			
Oxygen Scavengers	MEKO	19251	0.087 - 2.170 mg/l	MEKO	10+0	20	565	1.771			
Ozone	O ₃	00607	0.010 - 0.800 mg/l	O ₃	1+0	50	550	0.450			
Ozone	O ₃	00607	0.05 - 4.00 mg/l	O ₃	1+0	10	550	2.25			
Oxygen Demand, Biological	see BOD										
Oxygen Demand, Chemical	see COD										
Peroxide	see Hydrogen peroxide										

Programming data for Spectroquant[®] test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor
CT pH ¹⁾	pH	01744	6.4	-	8.8		pH	0+0	16	550	non linear
											factor a
											6.163
											factor b
											3.039
											factor c
											-1.913
											factor d
											0.393
											factor e
											0.082
Phenol	Phenol	00856	0.002	-	0.100	mg/l	Phenol	30+10	20	446	0.11
Phenol	Phenol	00856	0.025	-	1.000	mg/l	Phenol	10+0	50	495	1.67
Phenol	Phenol	00856	0.10	-	5.00	mg/l	Phenol	10+0	10	495	8.33
P(PMB)	P	14848	0.02	-	1.00	mg/l	P	5+0	50	690	0.417
P(PMB)	P	14848	0.10	-	5.00	mg/l	P	5+0	10	690	2.08
P(PMB)	PO4	14848	0.06	-	3.00	mg/l	PO4	5+0	50	690	1.275
P(PMB)	PO4	14848	0.3	-	15.3	mg/l	PO4	5+0	10	690	6.38
P(PMB)	P2O5	14848	0.05	-	2.30	mg/l	P2O5	5+0	50	690	0.952
P(PMB)	P2O5	14848	0.2	-	11.5	mg/l	P2O5	5+0	10	690	4.76
P(PMB)	PO4	14848	0.6	-	32.3	mmol/m3	PO4	5+0	50	690	13.4
P(PMB)	PO4	14848	3	-	161	mmol/m3	PO4	5+0	10	690	67.2
CT Phosphate	PO4-P	00616	3	-	100	mg/l	PO4-P	5+0	16	690	39.2
CT Phosphate	PO4	00616	9	-	307	mg/l	PO4	5+0	16	690	120
Phosphate	PO4-P	00798	1.0	-	100.0	mg/l	PO4-P	5+0	10	690	35.0
Phosphate	PO4	00798	3	-	307	mg/l	PO4	5+0	10	690	107

¹⁾ Remark to pH Cell Test:

As blank dist. water without adding reagents must be used.

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor
CT Potassium	K	14562	5.0	-	50.0	mg/l	K	5+0	16	690	non linear
										factor a	0
										factor b	50.586
										factor c	-24.882
										factor d	12.386
										factor e	0
CT Potassium	K	00615	30	-	300	mg/l	K	5+0	16	690	not linear
										factor a	0
										factor b	174.31
										factor c	-47.64
										factor d	23.35
										factor e	0
Silicate	Si	00857	0.5	-	50.0	mg/l	Si	2+2	10	405	25.3
Silicate	Si	00857	5	-	500	mg/l	Si	2+2	10	405	258
Silicate	SiO ₂	00857	1.1	-	107.0	mg/l	SiO ₂	2+2	10	405	54.3
Silicate	SiO ₂	00857	11	-	1070	mg/l	SiO ₂	2+2	10	405	550
Silver	Ag	14831	0.10	-	1.00	mg/l	Ag	5+0	50	550	0.444
Silver	Ag	14831	0.50	-	3.00	mg/l	Ag	5+0	10	550	2.22
CT Sodium in nutrient solutions for fertilizations	Na	00885	10	-	300	mg/l	Na	1+0	16	550	267

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor	
CT Sulfate	SO ₄	14548	5	-	250	mg/l	SO ₄	2+0	16	525	non linear	
											factor a	0
											factor b	132.55
											factor c	-75.27
											factor d	63.586
											factor e	0
CT Sulfate	SO ₄	00617	50	-	500	mg/l	SO ₄	2+0	16	525	non linear	
											factor a	0
											factor b	352.68
											factor c	-124.55
											factor d	137.57
											factor e	0
CT Sulfate	SO ₄	14564	100	-	1000	mg/l	SO ₄	2+0	16	525	non linear	
											factor a	0
											factor b	719.44
											factor c	-442.33
											factor d	369.43
											factor e	0
Sulfite	SO ₃	01746	1.0	-	60.0	mg/l	SO ₃	2+0	10	405	29.4	
Sulfite	SO ₃	01746	0.8	-	48.0	mg/l	SO ₃	2+0	10	405	23.5	
CT Surfactants (cationic)	CTAB	01764	0.05	-	1.50	mg/l	CTAB	5+0	16	635	1.92	
CT Surfactants (nonionic)	Triton X-100	01787	can not be measured in the SQ 118 because of not available wavelength									
Suspended solids	susS		25	-	750	mg/l	susS	0+0	20	820	730	
Tensides			see Surfactants									
Total Alkalinity			see Acid capacity to pH 4.3									

Programming data for Spectroquant[®] test kits for photometer SQ 118

Display-text		Item	Measuring range				Citation form	Reaction time	cells. mm	filter nm	factor
CT Total Hardness	Ca	00961	5	-	215	mg/l	Ca	3+0	16	565	non linear
										565	factor a
										565	factor b
										565	factor c
										565	factor d
										565	factor e
										565	0
										565	318.61
										565	-192.47
										565	76.5
										565	0
CT Total Hardness	CaCO ₃	00961	13	-	538	mg/l	CaCO ₃	3+0	16	565	non linear
										565	factor a
										565	factor b
										565	factor c
										565	factor d
										565	factor e
										565	0
										565	796.52
										565	-481.17
										565	191.29
										565	0
CT Total Hardness	CaO	00961	7	-	301	mg/l	CaO	3+0	16	565	non linear
										565	factor a
										565	factor b
										565	factor c
										565	factor d
										565	factor e
										565	0
										565	446.05
										565	-269.45
										565	107.12
										565	0
CT Total Hardness	°d	00961	0.7	-	30.1	mg/l	°d	3+0	16	565	non linear
										565	factor a
										565	factor b
										565	factor c
										565	factor d
										565	factor e
										565	0
										565	44.61
										565	-26.95
										565	10.71
										565	0

Programming data for Spectroquant® test kits for photometer SQ 118

Display-text	Item	Measuring range	Citation form	Reaction time	cells. mm	filter nm	factor	
CT Total Hardness °e	00961	0.9 - 37.6 mg/l	°e	3+0	16	565	non linear	
						factor a	0	
						factor b	55.76	
						factor c	-33.68	
						factor d	13.39	
						factor e	0	
CT Total Hardness °f	00961	1.3 - 53.8 mg/l	CaO	3+0	16	565	non linear	
						factor a	0	
						factor b	79.65	
						factor c	-48.12	
						factor d	19.13	
						factor e	0	
Total Nitrogen	see Nitrogen, total							
CT TOC	TOC	14878	5.0 - 80.0 mg/l	TOC	0+0	16	585	-64.5
CT TOC	TOC	14879	50 - 800 mg/l	TOC	0+0	16	585	-645
CT Volatile Organic Acids	HOAc	01763	50 - 3000 mg/l	HOAc	3+10	16	495	1430
Water hardness	see Residual hardness or Total hardness							
Zinc	Zn	14832	0.05 - 250 mg/l	Zn	3+2	10	565	1.08
CT Zinc	Zn	00861	0.025 - 1.000 mg/l	Zn	15+0	16	495	1.36

In our Internet homepage <http://photometry.merck.de> users can find all information about our test kits. All package inserts as well as manual pages from our NOVA or Pharo photometers can be found there. This information helps to get a quick overview about the easy handling of our test kits. The package inserts include detailed descriptions about all handling steps, colour reaction, sample preparation if required, interference with other substances, analytical quality assurance data and more.

Please contact our local Merck companies or Merck agents for any questions. Alternatively contact: environmental.analysis@merck.de.